



Butyl Persefose HP

1. Basic product information

The framework of Butyl Persefose HP hydrophobic resin is highly cross-linked agarose, which retains the excellent hydrophilicity and pore structure of natural polysaccharide compounds and has good compatibility with biological macromolecules. The average particle size of Butyl Persefose HP is 34 μm , which has a higher resolution efficiency than Butyl Persefose 4FF. The surface contains hydrophobic ligands; the ligand is butyl, which is particularly suitable for the fine separation and purification of recombinant proteins, antibodies, vaccines, virus-like particles and other highly hydrophobic biological molecules.

2. Chromatography resin parameters

Resin type	Hydrophobic interaction
Functional group	Butyl
Matrix	Highly cross-linked agarose
Median particle size	34 μm
Ligand density	~50 μmol butyl/ml
Recommended flow rate	60–120 cm/h
Maximum flow rate	200 cm/h
Maximum working pressure	3 bar
Working temperature	4–40°C

3. Chemical resistance

pH stability*	2–14 (CIP), 3–13 (working)
Chemical stability	All commonly used aqueous buffers, 30% isopropanol**, 70% ethanol**, 1M NaOH, 1M acetic acid, 6M guanidine hydrochloride

* The physical and chemical properties and functions of the chromatographic resin did not change significantly after being placed in an environment of 40°C and pH 3–13 for 7 days.

** v/v, volume ratio

4. Method of use

4.1 Chromatographic conditions

(1) The binding buffer is usually a phosphate buffer containing a high concentration of salt, such as 20–50 mM PB, 1.5–2 M $(\text{NH}_4)_2\text{SO}_4$, pH 7.0. The elution buffer normally uses phosphate buffer without other salts, such as 20–50 mM PB, pH 7.0, which needs to be based on the subsequent experimental results (whether there is precipitation of the target, binding strength of the target, recovery rate, resolution etc.) Adjust the concentration and type of salt in the binding buffer. For substances that are difficult to elute, pure water can be used, or low-concentration ethanol can be added to the pure water as the eluent.

(2) Flow rate: generally, choose a linear flow rate of 60–120 cm/h according to the column bed height.

(3) Sample pretreatment: to prevent the sample from clogging the column, the sample needs to be filtered with a 0.45 µm microporous membrane before loading it. It is recommended that the pH and conductivity of the sample are adjusted to be consistent with the equilibrium buffer (dilution, ultrafiltration can be used and desalting to adjust the pH and conductivity of the sample).

4.2 Chromatography steps

(1) Equilibration: use equilibration buffer to fully equilibrate the chromatography column until the pH and conductivity are stable and basically consistent with the equilibration buffer. This step usually requires 3–5 column bed volumes (CV).

(2) Sample loading*: according to the binding capacity measured in the small test, determine the sample loading volume and loading amount of the sample on Butyl Persefose HP.

(3) Impurity washing*: use equilibration buffer or other suitable buffer to wash the chromatography column until the UV stabilises and returns to the baseline.

(4) Elution*: elution is achieved by decreasing the concentration of salt ions. The concentration of salt ions in the elution buffer can be gradually decreased through a linear gradient or a step gradient to elute molecules with different binding strengths. pH gradient elution or mixed elution can also be used.

(5) Regeneration: rinse the column with a buffer containing low salt.

(6) Re-equilibration: re-equilibrate the column with equilibration buffer.

* Note: If the flow-through mode is used, the ‘sample loading’ step should be set to collect; the ‘washing’ step should ensure that the target molecules have completely flowed through, then the collection can be stopped; the ‘elution’ step should directly use low-salt buffer – impurities can be washed away.

5. Cleaning and regeneration

Contaminants (such as lipids, endotoxins and proteins) accumulate on the column as the number of uses of the chromatography resin increases. Determine the frequency of CIP according to the degree of contamination of the chromatography resin (if the contamination is considerable, CIP is recommended after each use to ensure repeatability of results and to prolong the working life of the chromatography resin). For different types of impurities and contaminants, the recommended cleaning conditions are as follows:


- Removal of strongly binding proteins: wash with 5 CV of ultrapure or pure water.
- Removal of strongly hydrophobic proteins and precipitated proteins: first wash with 5 CV of 1M NaOH solution, then wash the lye with 5–10 CV of ultrapure or pure water.
- Removal of lipoproteins and lipids: first wash with 5 CV of 70% ethanol or 30% isopropanol, then rinse with 5–10 CV of ultrapure or pure water.

Note: 70% ethanol or 30% isopropanol should be degassed before use; the flow rate should be 30–60 cm/h during CIP. Reverse cleaning should be used when the clogging is severe.

To reduce the microbial load, it is recommended that 0.5~1M NaOH solution is used to treat the chromatographic resin. Treatment time is 15~30 minutes.

6. Storage

Keep the unopened chromatography resin in the original container and store at 4~30°C in a well-ventilated, dry and clean place. Do not freeze. Wash the used column with 2–3 CV of 20% ethanol solution and store at 2~8°C.



LT BIOTECH LTD,
Mokslininkų g. 6A,
LT-08412 Vilnius,
Lithuania
www.ltbiotech.lt

7. Destruction and recycling

Since chromatography resin is difficult to degrade in nature, it is recommended that the waste chromatography resin is incinerated to protect the environment. For chromatography resin that has been in contact with biologically active samples such as viruses and blood, follow the local biosafety requirements before destroying or disposing of it.

8. Packing method

Detailed information on resin packaging is available on request. Please contact your local distributor.

9. Ordering information

Product name: Butyl Persefose HP

Product Cat. No	Package
136-00025	25 ml
136-00100	100 ml
136-00500	500 ml
136-01000	1 L
136-05000	5 L
136-10000	10 L
136-20000	20 L